

Unilateral Biportal Endoscopy Solution

BONSS RF Plasma Surgical System


UBE90



Hook-UBE



Spine-o-UBE



UXD70/UXD90

Unique angle and insulation design

Ensure safe operation, and optimized coagulation effect.
Designed for extraspinal operation.
Provide excellent clinical performance in dealing with different tissues.

Spine-o-UBE (135mm working length) Spine-o-QST (230mm working length)

Intraspinal use electrodes.
It is specifically designed for UBE when dealing with the tissue or hemostasis near the nerve root.
Its ergonomic hybrid handle allows for more control and usability with its variety of grip positions, as well as a wide custom comfort that conforms to your hand.



Based on our decades experience in R&D of RF Plasma technology, BONSS has developed the Unilateral Biportal Endoscopy (UBE) solution for minimally-invasive spine surgery. The patented electrode design is integrated with cutting, dissection, ablation, coagulation, hemostasis, and suction capabilities in one versatile device, and designed for intraspinal and extraspinal operations. BONSS RF Plasma Surgical System has become an integral part of UBE surgery.



UBE-AX

Dilator

Description: Dilator
Specification: Φ4mm(solid)*220mm
Model: UBE-A-KZGQ-2204

Features:

1. Establish both working channel and viewing channel
2. CT assisted positioning

Description: Semi-Tube
Specification: Φ10mm*50mm (L)
Model: UBE-A-KZGB-1060
Specification: Φ10mm*90mm (L)
Model: UBE-A-KZGB-1080
Specification: Φ12mm*60mm (L)
Model: UBE-A-KZGB-1010

Features:

1. Maintain working channel and smooth saline outlet
2. Ensure that the instrument enters the surgical site smoothly.

Description: Dilator
Specification: Φ10mm*120mm (L)
Model: UBE-A-KZGN-1611
Specification: Φ8mm*140mm (L)
Model: UBE-A-KZGN-1809
Specification: Φ6mm*160mm (L)
Model: UBE-A-KZGN-2007


Features:

1. Establishing the surgical channels and protect tissue.
2. Maintain surgical channels.
3. Expand the surgical field of view.



Instructions:

1. Insert the surgical site step by step after muscle stripping.
2. Place the tip on the surgical segment.



Instructions:

After the operation channel is established, the instrument tip is inserted into the body. The handpiece remains outside.



Instructions:

When using an instrument with a cutting edge, the instrument with a cutting edge is inserted into the body through the instrument channel to the operation site.

Tissue Retractor

Description: T-Retractor
Specification: Φ5mm(Thickness)*158mm(L)
Model: UBE-A-BLQA-2209
Features:For tissue detaching



Instructions:

After the surgical incision, insert its tip into the detaching position according to the scale of the instrument, swing it left and right, and strip the tissue.

Root Retractor

Description: Straight
Specification: length 100mm, width 4mm
Model: UBE-A-LGAL-1006
Specification: length 100mm, width 8mm
Model: UBE-A-LGAL-8006

Features:

1. For tissue traction.
2. For nerve traction.



Instructions:

After exposing the tissues that need to be pulled, extend the tip of the instrument and use the front hook to pull the tissues

Retractor

Description: ligamenta flava
Specification: length 100mm, widthh 2、3mm
Model: UBE-A-BLZB-0561
Specification: length 100mm, width 4、5mm
Model: UBE-A-BLZB-0562

Features:

It is used for ligamentum flavum traction, upgraded version.

Description: Double-ended
Specification: Edge width 3mm*angle 0°、15° *length 300mm
Model: UBE-A-BLZA-1504
Specification: Edge width 3mm*angle 25°、35° *length 300mm
Model: UBE-A-BLZA-2503

Features:

1. For tissue stripping
2. For anatomical exploration.



Instructions:

After exposing the tissues that need to be pulled, extend the tip of the instrument and use the front hook to pull the ligamentum flavum.



Instructions:

Separate two different tissues with the front flat part.

Protector

Description: Protector sheath
Specification: Φ9mm φ7mm*120mm (L) curved
Model: UBE-A-GLGW-9009

Features:

It is used to block tissue and expand the surgical field of view.



Instructions:

The instrument is inserted through the viewing channel with the tip against the end of the channel. Insert the assembled endoscope and sheath into the instrument. The front end shield can be used to shield tissue by rotating the instrument.

Guide Needle

Description: Flushing needle
Specification: φ4mm angle135°
Model: UBE-D-GDYZ-0016

Features:

Used to flush residual tissue in the intervertebral space.



Instructions:

The bottom is connected with the infusion tube, and the tip extends into the intervertebral space under the supervision of the endoscope, and the switch of the infusion tube is turned on for flushing.

Bone Probe

Description: L-shape
Specification: Length 100mm, Tip 4mm
Model: UBE-D-GTZA-1001

Features:

Used to explore nerves and confirm anatomical structures.



Instructions:

After the tissue is completely exposed, the tip of the instrument is inserted into the body for nerve exploration and confirmation of each tissue structure.

Retractor

Description: Disc Reamer
Specification: Length 220mm,Edge width 4mm
Model: UBE-A-BLZA-3503

Features:

Incision of the annulus.



Instructions:

After the annulus is exposed, the device is passed through the working channel to work on the annulus.

Pack box (UBE-AX)

Specification: (L)440mm*(W)290mm*(H)80mm



UBE-BX

Dilator

Description: Dilator
Specification: Φ12mm * 150mm(L)
Model: UBE-A-KZGS-1812

Specification: Φ14mm * 140mm(L)
Model: UBE-A-KZGS-1614

Specification: Φ16mm * 130mm(L)
Model: UBE-A-KZGS-2010

Features:
 Dilation of the surgical channel. For fusion surgery, it is necessary to expand the surgical channel.



Instructions:
 Insert the surgical site step by step after muscle stripping.

Chisel

Description: Upturned
Specification: Edge width 5mm (Upturned15°)
Model: UBE-B-GZZG-1204

Specification: Edge width 5mm (Upturned30°)
Model: UBE-B-GZQZ-1205

Description: Straight
Specification: Edge width 5mm, Straight
Model: UBE-B-GZQY-1205

Description: L-shape
Specification: Edge width 5mm, L-Type
Model: UBE-B-GZCZ-1205

Features:
 Chisel out the bone



Instructions:
 After the bony structure is exposed, put the cutting edge of the front end against the bone surface, and strike the tail of the osteotome with a bone hammer. The bevel faces the direction of bone resection.

Root Retractor

Description: Right bent
Specification: 100mm (L)*10mm (W)
Model: UBE-A-LGBY-8010

Description: Left bent
Specification: Edge width 5mm (Upturned30°)
Model: UBE-A-LGBZ-8010

Description: Right edgefold
Specification: 100mm (L)*5mm (W)
Model: UBE-A-LGRQ-1502

Description: Left edgefold
Specification: 100mm (L)*5mm (W)
Model: UBE-A-LGRQ-1501

Features:
 1. For tissue traction.
 2. For nerve traction.
 3. For CAGE implantation during fusion surgery.

Instructions:
 After exposing the tissues that need to be pulled, extend the tip of the instrument and use the front hook to pull the tissues



Retractor

Description: Cage Osteotome
Specification: Φ9mm(solid)*255mm
Model: UBE-A-BLZA-0505

Features:
 Used for CAGE position adjustment.



Instructions:
 After the CAGE is put in, insert the tip of the instrument into the threaded hole of the CAGE, and tap the bottom of the instrument with a bone hammer to adjust the position of the CAGE so that it lies horizontally.

Disc Reamer

Description: Straight
Specification: Edge width 7mm, length 200mm
Model: UBE-B-ZJDH-9505

Features:
 For reaming intervertebral discs.



Instructions:
 Insert the tip of the instrument into the intervertebral space, and rotate the handle to drive the cutting edge to rotate and ream the intervertebral disc.

Curette


Description: Straight
Specification: Edge width 3mm, length 140mm
Model: UBE-D-GCZA-1503

Specification: Edge width 4mm, length 140mm
Model: UBE-D-GCZA-1505

Description: Hollow
Specification: Edge width 3mm,length 140mm, Side-Bent
Model: UBE-D-GCZA-1506

Specification: Edge width 4mm,length 140mm, Side-Bent
Model: UBE-D-GCZA-1506

Features:
 1. For bone scraping
 2. Used to scrape endplates.



Instructions:
 Use the cutting edge of the working part of the instrument to scrape off the bone

Bone Hammer

Specification: 245mm(L)
Model: UBE-D-GCCZ-0450



Instructions:Used to knock other instruments.

Pack box (UBE-AX+BX)

Description: For UBE-AX and UBE-BX
Specification: (L)440mm*(W)290mm*(H)80mm



UBE-CX

Punch

Description: Dismounted rotatable

Specification: Edge width 2mm Length 220mm Angle 110°

Model: UBE-C-KTYB-1102

Specification: Edge width 3mm Length 220mm Angle 110°

Model: UBE-C-KTYB-1103

Specification: Edge width 2mm Length 220mm Angle 130°

Model: UBE-C-KTYB-1302

Specification: Edge width 3mm Length 220mm Angle 130°

Model: UBE-C-KTYB-1303

Description: Dismounted rotatable, curved

Specification: Edge width 2mm Length 220mm Angle 130°

Model: UBE-C-KTYA-1302

Specification: Edge width 3mm Length 220mm Angle 130°


Model: UBE-C-KTYA-1303

Features:

Rotatable, any angle can be adjusted

Instructions:

Biting and cutting the bone. Rotatable, any angle can be adjusted



Description: Handle for dismounted rotatable punch

Model: UBE-C-KTYB-0000

Features:

Handpiece

Instructions:

1. Align the arrow on the rongeur shaft with the handle arrow and insert into the handle;
2. Adjust the angle according to the needs of the operation, hold the bottom of the handle, press the rongeur rod in the direction of the arrow, rotate the rongeur rod, release the rongeur rod so that it is locked with the handle.
3. When disassembling, adjust the rongeur rod to align the arrow of the handle with the arrow of the rongeur rod according to the method in step 2, and then it can be disassembled.



Nucleus Pulposus Forceps

Description: Straight, Granular

Specification: Edge width 4mm Length 180mm

Model: UBE-C-SHZB-2004

Description: Upturned, Granular

Specification: Edge width 4mm Length 180mm

Model: UBE-C-SHWB-2004

Description: Upturned, Ball

Specification: Edge width 3mm Length 180mm

Model: UBE-C-SHWB-2003

Features:

1. Removal of tissue.
2. Used to remove intervertebral disc tissue.
3. The working part is reduced, which is more suitable for minimally invasive surgery.
4. The reduced diameter of the insertion section allows for easy access to and from the instrument channel.

Description: Straight, Square

Specification: Edge width 4mm Length 180mm

Model: UBE-C-SHZB-2003

Features:

Square bite forceps, neatly trim the ligamentum flavum and other soft tissues.



Pack box(UBE-CX)

Specification: 440mm (L)*290mm(W)*80mm(H)

UBE-EX

Dilator

Description: Scope sheath

Specification: Φ6mm * 120mm(L)* 0.5mm(Thickness), for 0° scope

Model: UBE-A-KZGJ-1206

Specification: Φ6mm * 115mm(L)* 0.5mm(Thickness) for 30° scope

Model: UBE-A-KZGJ-1506

Features:

1. Protect the endoscope.
2. Saline injection channel
3. Small and exquisite make is more ergonomic.



Bone Guide Needle

Description: Obturator

Specification: φ4mm * 130mm(L)

Model: UBE-D-GDYZ-1304

Model: UBE-D-GDYZ-1604

Features:

1. Protect the endoscope.
2. Saline injection channel
3. Small and exquisite make is more ergonomic.

Instructions:

1. After the endoscope is assembled with the sheath, it is inserted into the human body through the previously established observation channel.
2. The saline inlet is connected with the infusion tube for saline injection.

Pack box(E1X)

Pack box(E2X)

Description: One set of sheath in one box

Specification: 285mm(L)*100mm(W)*62mm(H)

