

Radiofrequency Ablation Portfolio

G4™ RF Generator



VERSATILE. INTUITIVE. RELIABLE.


The G4 RF Generator is designed to give consistent, reliable results with a wide range of features in an easy-to-use interface.

KEY FEATURES

- Intuitive 12" touchscreen with large color-coded status visualizations
- Ability to treat four sites at once with simultaneous or staggered start options
- Automatic or manual stimulation controlled with touchscreen or dial

TWO INTERFACES FOR CUSTOMIZED NEEDS

Each G4 user can create multiple customized presets.



ONE-TOUCH INTERFACE

- Large numeric readouts
- Simplified touchscreen controls




GRAPHING INTERFACE


- Graph temperature and voltage
- Dial or touchscreen control

FOUR INDEPENDENTLY CONTROLLED ELECTRODE OUTPUTS


Treat up to four sites at once using monopolar, bipolar, or pulsed RF for a wide range of treatment options.



KNEE JOINT




HIP JOINT




SACROILIAC JOINT




FACET JOINT




INTRADISCAL



INTRACRANIAL



TRIGEMINAL NEURALGIA




PODIATRY

DESIGNED FOR USE WITH BOSTON SCIENTIFIC RF PRODUCTS



ELECTRODES



CANNULAE



GROUNDING PADS

The Cosman Radiofrequency Generators, associated Radiofrequency Lesion Probes and RF Cannula are indicated for use in procedures to create radiofrequency lesions for the treatment of pain or for lesioning nerve tissue for functional neurosurgical procedures. The Cosman RF Injection Electrodes are used for percutaneous nerve blocks with local anesthetic solution or for radiofrequency lesioning. Refer to the Instructions for Use provided with Cosman generators, electrodes and cannulas for potential adverse effects, warnings and precautions prior to using these products.

All trademarks are the property of their respective owners.

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. Please check availability with your local sales representative or customer service.

Caution: Federal (U.S.) law restricts this device to sale by or on the order of a physician.

Outside of US Indications for Use: CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for use only in countries with applicable health authority registrations. Material not intended for use in France.

**Boston
Scientific**
Advancing science for life™

25155 Rye Canyon Loop
Valencia, CA 91355 USA

Copyright ©2017 by
Boston Scientific Corporation
or its affiliates. All rights reserved.

NM-328703-AA_JAN2017